



27329

## o-Phthalaldehyde (OPA) extrapure AR, 99%

Part B

CAS : 643-79-8

Molecular Formula : C<sub>8</sub>H<sub>6</sub>O<sub>2</sub>

Molecular Weight : 134.13

### Specifications

Appearance (Colour)	Pale yellow
Appearance (Form)	Crystalline powder
Solubility (Turbidity) 10% solution in methanol	Clear
Solubility (Colour) 10% solution in methanol	Yellow
Assay	min. 99%
Melting Point	54 - 56°C
Sulphated Ash	max. 0.1%
Acidity	max. 1% (Phthalic acid)

### Other Information

Application

Reagent for amino acids Used for GC Derivatization.

### General Information

Storage	2 to 8°C (Refrigerate)
Shelf Life	60 Months
IMDG Identification	UN No.:2923 , IMCO Class No.:8(6.1) , Packing Group:II
HSN Code	
25 Gms	29122990 (GST 18%)
5 Gms	29122990 (GST 18%)
Type of Packing	
25 Gms	Glass Bottle

### Available Packages

5 Gms

25 Gms

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com