



27391

Chitosan (Low MW) extrapure, 10-150m.Pas, 90% DA

Part B

CAS : 9012-76-4

Specifications

Appearance (Colour)	White to light yellow
Appearance (Form)	Free flowing powder
Solubility (Colour) 1% solution in 1% Acetic acid	Colourless to pale yellow
Loss on drying	max. 10%
Viscosity (20°C)	10 - 150 mPa.s
Deacetylated Degree	min. 90%
Not for medicinal use	.

Other Information

Application

Biocompatible, antibacterial and environmentally friendly polyelectrolyte with a variety of applications including water treatment, chromatography, additives for cosmetics, textile treatment for antimicrobial activity, novel fibers for textiles, photographic papers, biodegradable films, biomedical devices, and micro-capsule implants for controlled release in drug delivery.

General Information

Storage	8 to 25°C (Cool & Dry Area)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
100 Gms	2921190 (GST 18%)
25 Gms	2921190 (GST 18%)
500 Gms	2921190 (GST 18%)
Type of Packing	
100 Gms	Glass Bottle
25 Gms	Glass Bottle
500 Gms	Plastic Bottle

Available Packages

25 Gms

100 Gms

500 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com
