



27587

Glycyl Leucine Phenylalanine extrapure, 99%

Part B

CAS : 103213-38-3

Molecular Formula : C₁₇H₂₅N₃O₄

Molecular Weight : .335.40

Other Information

Application

Immunostimulating peptide from casein.

General Information

Storage

2 to 8°C (Refrigerate)

Shelf Life

60 Months

IMDG Identification

Not Regulated for Transport (Non-Haz)

HSN Code

100 Mg

29379020 (GST 18%)

Available Packages

100 Mg

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com