



27821

## Silver Trifluoromethanesulfonate extrapure, 98%

Part B

CAS : 2923-28-6

Molecular Formula : CF<sub>3</sub>SO<sub>3</sub>Ag

Molecular Weight : 256.94

### Specifications

Appearance (Colour)	White to off-white or greyish
Appearance (Form)	compound
Assay	min. 98%

### General Information

Storage	2 to 8°C (Refrigerate)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
1 Gms	28432900 (GST 18%)
25 Gms	28432900 (GST 18%)
5 Gms	28432900 (GST 18%)

### Available Packages

1 Gms

5 Gms

25 Gms

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com