



27956

## Dithioerythritol (DTE) extrapure, 99%

Part B

CAS : 6892-68-8

Molecular Formula : C<sub>4</sub>H<sub>10</sub>O<sub>2</sub>S<sub>2</sub>

Molecular Weight : 154.25

### Specifications

Appearance (Colour)	White
Appearance (Form)	Crystalline Compound
Solubility (Turbidity) 5% aq. solution	Clear
Solubility (Colour) 5% aq. solution	Colourless
Assay	min. 99%
Melting Point	82 - 84°C

### Other Information

Application	Reagent for cleavage of disulphide bonds in proteins.
-------------	---

### General Information

Storage	0 to 4°C (Freeze)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
1 Gms	29309099 (GST 18%)
25 Gms	29309099 (GST 18%)
5 Gms	29309099 (GST 18%)
Type of Packing	
5 Gms	Glass Bottle

### Available Packages

1 Gms

5 Gms

25 Gms

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com