



27999

Asiatic Acid pure, 97%

Part B

CAS : 464-92-6

Molecular Formula : C₃₀H₄₈O₅

Molecular Weight : 488.70

Specifications

Appearance (Colour)	White to off - white
Appearance (Form)	Powder
Assay (HPLC)	min. 97%

Other Information

Application

Asiatic acid is commonly used in wound healing. Asiatic acid has antioxidant, anti-inflammatory and neuroprotective properties.

General Information

Storage	0 to 4°C (Freeze)
Shelf Life	60 Months
HSN Code	
200 Mg	29161990 (GST 18%)
50 Mg	29161990 (GST 18%)

Available Packages

50 Mg
200 Mg

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com