



28033

L-Thioproline extrapure, 98%

Part B

CAS : 34592-47-7

Molecular Formula : C₄H₇NO₂S

Molecular Weight : 133.17

Specifications

Appearance (Colour)	White to off - white
Appearance (Form)	Crystalline powder
Assay	min. 98%

Other Information

Application	Used in peptide coupling reactions.
-------------	-------------------------------------

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
HSN Code	
10 Gms	29341000 (GST 18%)
25 Gms	29341000 (GST 18%)

Available Packages

10 Gms
25 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com