



28104

Phenol Crystalline pure, 99%

Part A

CAS : 108-95-2

Molecular Formula : C₆H₆O

Molecular Weight : 94.11

Specifications

Appearance (Colour)	White to slight yellow
Appearance (Form)	Crystalline solid
Assay (GC)	min. 99%
Melting Point	40-41°C
Water (KF)	max. 0.5%

General Information

Storage	8 to 25°C (Cool & Dry Area)
Shelf Life	60 Months
IMDG Identification	UN No.:1671 , IMCO Class No.:6.1 , Packing Group:II
HSN Code	
2.5 Kg	29071110 (GST 18%)
25 Kg	29071110 (GST 18%)
500 Gms	29071110 (GST 18%)
Type of Packing	
2.5 Kg	Glass Bottle

Available Packages

500 Gms

2.5 Kg

25 Kg

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com