



28354

Dextran Blue 5 ex. Leuconostoc Sp.

Part B

CAS : 87915-38-6

Specifications

Appearance (Colour)	Blue to dark blue
Appearance (Form)	Powder
Average Molecular weight	made from Dextran 5 KDa
Loss on drying	max. 10%
Degree of substitution	0.01 - 0.08

Other Information

Application

Blue dextran is used in affinity chromatography, gel filtration chromatography, protein chromatography and molecular weight markers. Blue dextran has been used to study Alzheimer's disease.

General Information

Storage	-20 °C (Blue/Dry Ice)
Shelf Life	60 Months
HSN Code	
500 Mg	35051090 (GST 18%)

Available Packages

500 Mg

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com