



28599

Polysorbate 20 (Tween 20) extrapure

Part A

CAS : 9005-64-5

Specifications

Appearance (Clarity)	Clear
Appearance (Colour)	Pale yellow
Appearance (Form)	Viscous liquid
Density (g/ml) @ 20°C	1.105 - 1.107 g/ml
Refractive Index (20°C)	1.470 - 1.472
Acid value	max. 2.0
Saponification value	40 - 50
Hydroxyl value	96 - 108
Water (KF)	max. 2%
IR Identification Test	Passes Test

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
2.5 Kg	34024900 (GST 18%)
25 Kg	34024900 (GST 18%)
5 Kg	34024900 (GST 18%)
500 Gms	34024900 (GST 18%)
Type of Packing	
500 Gms	Glass Bottle

Available Packages

500 Gms
2.5 Kg
5 Kg
25 Kg

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com