



28682

D-(-)-Tartaric Acid extrapure, 99%

Part B

CAS : 147-71-7

Molecular Formula : C₄H₆O₆

Molecular Weight : 150.09

Specifications

Appearance (Colour)	White to off - white
Appearance (Form)	Crystalline powder
Solubility (Turbidity) 10% aq. solution	Clear
Solubility (Colour) 10% aq. solution	Colourless
Assay (by Neutralization)	min. 99%
Specific rotation [?]	-13 to -14° (C=10,water)
Melting Point	166 - 169°C
Sulphated Ash	max. 0.05%

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
100 Gms	29181200 (GST 18%)
250 Gms	29181200 (GST 18%)
Type of Packing	
100 Gms	Plastic Bottle

Available Packages

100 Gms

250 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com