



28699

## Selenium Nanopowder in aq. media, 0.1mg/ml, 50nm

Part C

CAS : 7782-49-2

Atomic Formula : Se

Atomic Weight : 78.96

### Specifications

Appearance (Colour)	Red
Appearance (Form)	Suspension
Concentration	0.1 mg/ml
APS	50 nm

### General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:3077 , IMCO Class No.:9 , Packing Group:III
HSN Code	
2 ml	28049000 (GST 18%)

### Available Packages

2 ml

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com