



29444

## Bismuth AAS Standard Solution in 0.5N HNO<sub>3</sub> (Traceable to NIST)

Part A

### Specifications

Appearance (Clarity)	Clear
Appearance (Colour)	Colourless
Appearance (Form)	Liquid
Concentration (Bi)	1000 ± 50 mg/Ltr
Traceability	Traceable to NIST

### General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	36 Months
IMDG Identification	UN No.:3264 , IMCO Class No.:8 , Packing Group:III
HSN Code	
100 ml	38229090 (GST 12%)
250 ml	38229090 (GST 12%)

### Available Packages

100 ml

---

250 ml

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com