



29467

Potassium Hydrogen Phthalate extrapure AR, ACS, ExiPlus, Multi-Compendial, 99.9%

Part A

CAS : 877-24-7

Molecular Formula : C₈H₅KO₄

Molecular Weight : 204.23

Specifications

ExiPlus	meets compendial specs of BP, Ph.Eur
Appearance (Colour)	White
Appearance (Form)	Crystalline powder
Solubility (Turbidity) 10% aq. solution	Clear
Solubility (Colour) 10% aq. solution	Colourless
Assay (by Neutralization)	min. 99.9%
pH (5% aq. solution)	4.0 - 4.02
Insoluble matter	max. 0.005%
Loss on drying	max. 0.05%
Chloride (Cl)	max. 0.002%
Sulphate (SO ₄)	max. 0.005%
Iron (Fe)	max. 0.0005%
Heavy Metals (Pb)	max. 0.0005%
Sulphur Compound (S)	max. 0.002%
Phthalic acid	max. 0.005%
Sodium (Na)	max. 0.005%

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
500 Gms	29173990 (GST 18%)
Type of Packing	
500 Gms	Plastic Bottle

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com