



29563

## Strontium Nitrate Anhydrous extrapure, 98%

Part A

CAS : 10042-76-9

Molecular Formula :  $\text{Sr}(\text{NO}_3)_2$

Molecular Weight : 211.63

### Specifications

Appearance (Colour)	White
Appearance (Form)	Crystalline compound / powder
Solubility (Turbidity) 10% aq. solution	Clear
Solubility (Colour) 10% aq. solution	Colourless
pH (5% aq. solution)	5.0 - 7.0
Assay (EDTA Titration)	min. 98%
Chloride (Cl)	max. 0.005%
Sulphate (SO <sub>4</sub> )	max. 0.01%
Iron (Fe)	max. 0.001%
Lead (Pb)	max. 0.001%

### General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:1507 , IMCO Class No.:5.1 , Packing Group:III
HSN Code	
500 Gms	28500049 (GST 18%)
Type of Packing	
500 Gms	Plastic Bottle

### Available Packages

500 Gms

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com