



29740

Tetramethylammonium Hydroxide Pentahydrate extrapure, 98%

Part A

CAS : 10424-65-4

Molecular Formula : (CH₃)₄N(OH)·5H₂O

Molecular Weight : 181.23

Specifications

| | |
|-----------------------------|----------------------|
| Appearance (Colour) | White |
| Appearance (Form) | Crystalline compound |
| Assay (by Neutralization) | min. 98% |
| Chloride (Cl) | max. 0.1% |
| Sulphate (SO ₄) | max. 0.1% |

Other Information

Notes

Used in Quality Control Testing for
Pharmaceutical API's and
Formulations.

General Information

| | |
|---------------------|--|
| Storage | 25 to 40°C (Room Temperature) |
| Shelf Life | 60 Months |
| IMDG Identification | UN No.:3423 , IMCO Class No.:8 , Packing Group:II |
| HSN Code | |
| 100 Gms | 29239000 (GST 18%) |
| 25 Gms | 29239000 (GST 18%) |
| 5 Gms | 29239000 (GST 18%) |
| Type of Packing | |
| 100 Gms | Glass Bottle |
| 25 Gms | Glass Bottle |

Available Packages

5 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com