



29821

Sodium Fluoride extrapure AR, 99%

Part A

CAS : 7681-49-4

Molecular Formula : NaF

Molecular Weight : 41.99

Specifications

Appearance (Colour)	White
Appearance (Form)	Crystalline powder
Solubility (Turbidity) 2.5% hot aq. solution	Clear
Solubility (Colour) 2.5% hot aq. solution	Colourless
Assay	min. 99%
Loss on drying	max. 0.2%
Acidity (HF)	max. 0.05%
Sulphate (SO ₄)	max. 0.01%
Iron (Fe)	max. 0.002%
Heavy Metals (Pb)	max. 0.001%
Copper (Cu)	max. 0.0005%
Sulphite (SO ₃)	max. 0.005%
Fluorosilicate (SiF ₆)	max. 0.1%

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Restricted for export from India
HSN Code	
5 Kg	28261990 (GST 18%)
500 Gms	28261990 (GST 18%)
Type of Packing	
500 Gms	Plastic Bottle

Available Packages

500 Gms

5 Kg

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com