



30064

Tetrabutylammonium Hydroxide 54-56% aq. solution extrapure AR

Part A

CAS : 2052-49-5

Molecular Formula : C₁₆H₃₇NO

Molecular Weight : 259.48

Specifications

Appearance (Clarity)	Clear
Appearance (Colour)	Colourless to pale yellow
Appearance (Form)	Liquid
Assay (T)	54-56%
Density (g/ml) @ 20°C	0.968-0.969
Halides (Br)	max. 0.1%
Sulphate (SO ₄)	max. 0.2%

Other Information

Notes

Highly concentrated Quaternary Ammonium solution, used in Quality Control Testing for Pharmaceutical API's and Formulations.

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:3267 , IMCO Class No.:8 , Packing Group:II
HSN Code	
100 ml	38229090 (GST 12%)
500 ml	38229090 (GST 12%)
Type of Packing	
500 ml	Glass Bottle

Available Packages

100 ml

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com