



30163

L-Lysine (free base) Anhydrous extrapure, 98%

Part A

CAS : 56-87-1

Molecular Formula : C₆H₁₄N₂O₂

Molecular Weight : 146.20

Specifications

Appearance (Colour)	White to pale yellow
Appearance (Form)	Crystalline powder
Solubility (Turbidity) 10% solution 2N HCl	Clear
Solubility (Colour) 10% solution 2N HCl	Colourless to pale yellow
Assay (NT)	min. 98%
Melting Point	213 - 214°C
Specific rotation [?]	+25° to +26° [C=5 , 5N HCL]
Ammonium (NH ₄)	max. 0.01%
Water (KF)	max. 2%

General Information

Storage	2 to 8°C (Refrigerate)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
10 Gms	29224100 (GST 18%)
100 Gms	29224100 (GST 18%)
25 Gms	29224100 (GST 18%)
Type of Packing	
25 Gms	Glass Bottle

Available Packages

10 Gms
25 Gms
100 Gms
500 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com