



30315

Tetrabutylammonium Fluoride (TBAF) 1M Soln. in THF

Part A

CAS : 429-41-4

Molecular Formula : C₁₆H₃₆FN

Molecular Weight : 261.47

Specifications

Appearance (Clarity)	Clear
Appearance (Colour)	Pale yellow to brown
Appearance (Form)	Liquid
Molarity (Titration by HClO ₄)	0.95 - 1.10 M
Density (g/ml) @ 25°C	0.907 - 0.909
Water (KF)	4.5 - 5.0%

General Information

Storage	2 to 8°C (Refrigerate)
Shelf Life	24 Months
IMDG Identification	UN No.:1993 , IMCO Class No.:3 , Packing Group:II
HSN Code	
100 ml	38220090 (GST 12%)
500 ml	38220090 (GST 12%)

Available Packages

100 ml

500 ml

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srlchem.com, website www.srlchem.com