



30565

Xylan from Corn Cob, 95%

Part B

CAS : 9014-63-5

Specifications	
Appearance (Colour)	White to off-white
Appearance (Form)	Crystalline powder
Assay	min. 95%
pH	3.5 - 6.5
Loss on drying	max. 5.0%

General Information	
Storage	8 to 25°C (Cool & Dry Area)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
10 Gms	29400000 (GST 18%)
100 Gms	29400000 (GST 18%)
25 Gms	29400000 (GST 18%)
500 Gms	29400000 (GST 18%)
Type of Packing	
10 Gms	Plastic Bottle
100 Gms	Plastic Bottle
500 Gms	Plastic Bottle

Available Packages	
10 Gms	
25 Gms	
100 Gms	
500 Gms	

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srlchem.com