



30795

Nickel (II) Nitrate Hexahydrate extrapure AR, 99%

Part A

CAS : 13478-00-7

Molecular Formula : Ni(NO₃)₂.6H₂O

Molecular Weight : 290.79

Specifications

Appearance (Colour)	Green
Appearance (Form)	Crystalline powder / Deliquescent crystals
Solubility (Turbidity) 10% aq. solution	Clear
Assay (EDTA Titration)	min. 99%
Chloride (Cl)	max. 0.001%
Sulphate (SO ₄)	max. 0.005%
Iron (Fe)	max. 0.001%
Ammonium (NH ₄)	max. 0.05%

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:2725 , IMCO Class No.:5.1 , Packing Group:III
HSN Code	
250 Gms	28342990 (GST 18%)
500 Gms	28342990 (GST 18%)
Type of Packing	
250 Gms	Glass Bottle
500 Gms	Glass Bottle

Available Packages

250 Gms
500 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com