



30888

**Sodium Tetrphenyl Borate (Kalignost)**  
**extrapure AR, ACS, ExiPlus, Multi-Compendial,**  
**99.5%**

Part A

CAS : 143-66-8

Molecular Formula : C<sub>24</sub>H<sub>20</sub>BNa

Molecular Weight : 342.23

### Specifications

ExiPlus	meets compendial specs of BP, Ph.Eur
Appearance (Colour)	White to beige
Appearance (Form)	Solid
Assay	min. 99.5%
Loss on drying	max. 0.5%
Chloride (Cl)	max. 0.005%
Sulphate (SO <sub>4</sub> )	max. 0.005%
Heavy Metals (Pb)	max. 0.001%
Clarity of solution	Passes test

### General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:2811 , IMCO Class No.:6.1 , Packing Group:III
HSN Code	
10 Gms	29319090 (GST 18%)
100 Gms	29319090 (GST 18%)
Type of Packing	
100 Gms	Glass Bottle

### Available Packages

10 Gms  
100 Gms  
1 Kg

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com