



31018

Cedarwood Oil Natural for microscopy

Part A

CAS : 8000-27-9

Specifications	
Appearance (Clarity)	Clear
Appearance (Colour)	Colourless to yellow to orange to amber
Appearance (Form)	Oily liquid
Density (g/ml) @ 20°C	0.930 - 0.940 g/ml
Refractive Index (20°C)	1.513 - 1.519

Other Information	
Notes	We also offer Immersion Oil Synthetic grade.

General Information	
Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
125 ml	33012916 (GST 18%)
30 ml	33012916 (GST 18%)
500 ml	33012916 (GST 18%)

Available Packages	
30 ml	
125 ml	
500 ml	

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srlchem.com