



31196

## Magnesium Chloride Anhydrous extrapure, 98%

Part A

CAS : 7786-30-3

Molecular Formula : MgCl<sub>2</sub>

Molecular Weight : 95.21

### Specifications

Appearance (Colour)	White
Appearance (Form)	Deliquescent flakes ,may contain some dark particles
Assay (EDTA Titration)	min. 98%
Loss on drying	max. 2%
Sulphate (SO <sub>4</sub> )	max. 0.02%
Heavy Metals (Pb)	max. 0.002%

### General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
100 Gms	28273100 (GST 18%)
5 Kg	28273100 (GST 18%)
500 Gms	28273100 (GST 18%)
Type of Packing	
100 Gms	Plastic Bottle
5 Kg	Plastic Drum
500 Gms	Plastic Bottle

### Available Packages

100 Gms  
500 Gms  
5 Kg

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com