



31398

a-Benzoinoxime extrapure AR, 99%

Part B

CAS : 441-38-3

Molecular Formula : C₁₄H₁₃NO₂

Molecular Weight : 227.26

Specifications

Appearance (Colour)	White to off - white
Appearance (Form)	Crystalline powder
Solubility (Turbidity) 2% solution in Methanol	Clear
Solubility (Colour) 2% solution in Methanol	Almost colourless to pale yellow
Assay (ex N)	min. 99%
Melting Point	153 - 155°C
Loss on drying	max. 0.5%
Sulphated Ash	max. 0.05%
Sensitivity to copper (Cu)	Passes test
Sensitivity to Molybdenium	Passes test

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
100 Gms	29420090 (GST 18%)
25 Gms	29420090 (GST 18%)
Type of Packing	
100 Gms	Plastic Bottle
25 Gms	Plastic Bottle

Available Packages

25 Gms
100 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com