



31499

2,4-Dichlorophenoxyacetic Acid (2,4-D) for tissue culture, 98%

Part A

CAS : 94-75-7

Molecular Formula : C₈H₆Cl₂O₃

Molecular Weight : 221.04

Specifications

Appearance (Colour)	White to off-white
Appearance (Form)	Powder
Solubility (Turbidity) 5% solution in Methanol	Clear
Assay (NaOH titration)	min. 98%
Melting Point	136 - 140°C
Suitability for Tissue Culture	Passes

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Restricted for export from India
HSN Code	
100 Gms	98020000 (GST 18%)
25 Gms	98020000 (GST 18%)
500 Gms	98020000 (GST 18%)
Type of Packing	
100 Gms	Glass Bottle
500 Gms	Glass Bottle

Available Packages

25 Gms
100 Gms
500 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com