



31751

Potassium Chromate extrapure AR, ACS, ExiPlus, Multi-Compendial, 99.5%

Part A

CAS : 7789-00-6

Molecular Formula : K₂CrO₄

Molecular Weight : 194.19

Specifications

ExiPlus	meets compendial specs of USP, BP, Ph.Eur
Appearance (Colour)	Yellow
Appearance (Form)	Crystalline powder
Solubility (Turbidity) 10% aq. solution	Clear
Assay (Iodometric)	min. 99.5%
pH (5% aq. solution)	9.0 - 9.8
Insoluble matter	max. 0.005%
Chloride (Cl)	max. 0.001%
Sulphate (SO ₄)	max. 0.01%
Sodium (Na)	max. 0.02%
Calcium (Ca)	max. 0.005%

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:3288 , IMCO Class No.:6.1 , Packing Group:III
HSN Code	
500 Gms	28415090 (GST 18%)
Type of Packing	
500 Gms	Plastic Bottle

Available Packages

500 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com