



31904

## Safranin O for cell culture, Endotoxin (BET) 0.05EU/mg

Part B1

CAS : 477-73-6

Molecular Formula : C<sub>20</sub>H<sub>19</sub>N<sub>4</sub>Cl

Molecular Weight : 350.84

### Specifications

Appearance (Colour)	Dark green
Appearance (Form)	Crystalline powder
Absorbance (A) of 1% aq. solution in a 1cm cell v/s H <sub>2</sub> O	
@520nm	min. 900
Loss on drying	max. 5%
Suitability for Cell Culture	Passes test
Bacterial Endotoxin Test (BET) (EU/mg = IU/mg)	0.05 EU/mg

### General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
100 Gms	32049000 (GST 18%)
25 Gms	32049000 (GST 18%)
Type of Packing	
100 Gms	Glass Bottle
25 Gms	Glass Bottle

### Available Packages

25 Gms  
100 Gms

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srlchem.com, website www.srlchem.com