



32212

Peppermint oil extrapure, 50%

Part B

CAS : 8006-90-4

Specifications

Appearance (Clarity)	Clear
Appearance (Colour)	Slight pale yellow
Appearance (Form)	Liquid
Assay (GC)	min. 50%

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
100 ml	33012400 (GST 18%)
500 ml	33012400 (GST 18%)

Available Packages

100 ml
500 ml

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com