



32431

tert-Butylchlorodiphenylsilane (TBDPSCI)
extrapure AR, 98%

Part B

CAS : 58479-61-1

Molecular Formula : C16H19ClSi

Molecular Weight : 274.86

Specifications

Appearance (Clarity)	Clear
Appearance (Colour)	Colourless to slight yellow
Appearance (Form)	Liquid
Assay (GC)	min. 98%
Density (g/ml) @ 20°C	1.072-1.074
Refractive Index (20°C)	1.566-1.568

General Information

Storage	2 to 8°C (Refrigerate)
Shelf Life	60 Months
IMDG Identification	UN No.:2987 , IMCO Class No.:8 , Packing Group:II
HSN Code	
10 ml	29319090 (GST 18%)
100 ml	29319090 (GST 18%)
25 ml	29319090 (GST 18%)
Type of Packing	
25 ml	Glass Bottle

Available Packages

10 ml
25 ml
100 ml

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srlchem.com