



32701

## tert-Butyl Hydroquinone (TBHQ) extrapure AR, 98%

Part A

CAS : 1948-33-0

Molecular Formula : C<sub>10</sub>H<sub>14</sub>O<sub>2</sub>

Molecular Weight : 166.22

### Specifications

Appearance (Colour)	White to light brown
Appearance (Form)	Powder/ crystals
Solubility (Turbidity) 10% solution in methanol	Clear
Assay (GC)	min. 98%
Melting Point	125 - 131°C

### General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:3077 , IMCO Class No.:9 , Packing Group:III
HSN Code	
2.5 Kg	29147990 (GST 18%)
250 Gms	29147990 (GST 18%)
Type of Packing	
250 Gms	Plastic Bottle

### Available Packages

250 Gms  
2.5 Kg

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com