



32737

## Silver AAS Standard Solution in 0.5N HNO<sub>3</sub> (Traceable to NIST)

Part A

### Specifications

|                      |                   |
|----------------------|-------------------|
| Appearance (Clarity) | Clear             |
| Appearance (Colour)  | Colourless        |
| Appearance (Form)    | Liquid            |
| Concentration (Ag)   | 1000 ± 50 mg/Ltr  |
| Traceability         | Traceable to NIST |

### General Information

|                     |   |
|---------------------|---|
| Storage             | 25 to 40°C (Room Temperature)                         |
| Shelf Life          | 36 Months   |
| IMDG Identification | UN No.:3264 , IMCO Class No.:8 ,<br>Packing Group:III |
| HSN Code            |   |
| 100 ml              | 38229090 (GST 12%)                                    |
| 250 ml              | 38229090 (GST 12%)                                    |

### Available Packages

100 ml

---

250 ml

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com