



32900

Z-DL-Phenylalanine extrapure, 99%

Part B

CAS : 3588-57-6

Molecular Formula : C₁₇H₁₇NO₄

Molecular Weight : 299.30

General Information

Storage	0 to 4°C (Freeze)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
25 Gms	29379020 (GST 18%)
5 Gms	29379020 (GST 18%)
Type of Packing	
25 Gms	Plastic Bottle

Available Packages

5 Gms

25 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com