



32915

Bismuth (III) Nitrate Basic extrapure AR, 79% BiNO₃

Part A

CAS : 1304-85-4

Molecular Formula : H₉Bi₅N₄O₁₂

Molecular Weight : 1461.99

Specifications

Appearance (Colour)	White to almost white
Appearance (Form)	Powder
Assay (Bi ₂ O ₃)	min. 79%
Content	70 to 74% (Bi)
Loss on drying	max. 3%

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:1477 , IMCO Class No.:5.1 , Packing Group:II
HSN Code	
100 Gms	28342990 (GST 18%)
500 Gms	28342990 (GST 18%)
Type of Packing	
100 Gms	Glass Bottle
500 Gms	Glass Bottle

Available Packages

100 Gms
500 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srlchem.com, website www.srlchem.com