



33138

Exton`s Reagent

Part B

Specifications

Appearance (Clarity)	Clear
Appearance (Colour)	Faint pink
Appearance (Form)	Solution
Suitability as reagent for qualitative detection of albumin	Passes test

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	36 Months
HSN Code	
100 ml	38229090 (GST 12%)

Available Packages

100 ml

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com