



33203

## Mid-Range 1 Protein Marker: 14-80 kDa

Part E

### Specifications

Appearance (Colour)	Pure blue
Appearance (Form)	Liquid
Appearance (Clarity)	Clear
SDS-PAGE Analysis	Electrophoresis of the marker should show nine crisp bands of almost equal intensity and be visualised by Coomassie Blue staining or SpryBlue staining.

### Other Information

Description	BioLit Midrange 1 Protein Marker covers the useful range of protein sizes that a researcher or student requires for routine experiments. This is a ready-to-use protein marker. The markers have been carefully prepared under sterile conditions and stored in appropriate additives to ensure long term stability.
Application	SDS PAGE Analysis & Western Blotting

### General Information

Storage	-20 °C (Blue/Dry Ice)
Shelf Life	18 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
0.5 ml	38229090 (GST 12%)
1 ml	38229090 (GST 12%)
2.5 ml	38229090 (GST 12%)

### Available Packages

0.5 ml

1 ml

2.5 ml

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com