



33340

Cellulose microcrystalline pure for TLC, 98%

Part B

CAS : 9004-34-6

Specifications

Appearance (Colour)	White to off-white
Appearance (Form)	Powder
pH (2% aq. solution)	5.0 - 7.5
Assay	min. 98%
Loss on drying	max. 6%

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
25 Kg	39129090 (GST 18%)
5 Kg	39129090 (GST 18%)
500 Gms	39129090 (GST 18%)
Type of Packing	
500 Gms	Plastic Bottle

Available Packages

500 Gms

5 Kg

25 Kg

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com