



33401

## D-Isoleucine extrapure, 98%

Part B

CAS : 319-78-8

Molecular Formula : C<sub>6</sub>H<sub>13</sub>NO<sub>2</sub>

Molecular Weight : 131.17

### Specifications

Appearance (Colour)	White to off - White
Appearance (Form)	Crystalline powder
Assay	min. 98%

### General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
1 Gms	29224990 (GST 18%)
250 Mg	29224990 (GST 18%)
5 Gms	29224990 (GST 18%)

### Available Packages

250 Mg

1 Gms

5 Gms

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com