



33434

DEAE-Dextran Hydrochloride 500 ex.
Leuconostoc Sp.

Part B

CAS : 9064-91-9

Specifications

Appearance (Colour)	White to off-white
Appearance (Form)	Powder
Loss on drying	max. 10%
Nitrogen (N) content	3.0 - 3.8%
pH (5% aq. solution)	5.0 - 7.0
Sulphated Ash	max. 2%
Average Molecular weight	made from Dextran 400 - 750 KDa
Free DEAE-Cl	max. 0.002%

General Information

Storage	-20 °C (Blue/Dry Ice)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
100 Gms	39140090 (GST 18%)
25 Gms	39140090 (GST 18%)
5 Gms	39140090 (GST 18%)

Available Packages

5 Gms
25 Gms
100 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srlchem.com