



34103

## DPPP (1,3-Bis(Diphenylphosphino)Propane) extrapure, 98%

Part B

CAS : 6737-42-4

Molecular Formula : C<sub>27</sub>H<sub>26</sub>P<sub>2</sub>

Molecular Weight : 412.44

### Specifications

Appearance (Colour)	White to off-white
Appearance (Form)	Solid
Assay	min. 98%
Melting Point	61 - 64 °C

### Other Information

Application

Generally used in the following chemical reactions namely C-C Bond Formation, Carbonylations, Negishi Coupling, Sonogashira Coupling, Suzuki-Miyaura Coupling.

### General Information

Storage	8 to 25°C (Cool & Dry Area)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
1 Gms	29029090 (GST 18%)
25 Gms	29029090 (GST 18%)
5 Gms	29029090 (GST 18%)
Type of Packing	
1 Gms	Plastic Bottle
25 Gms	Glass Bottle
5 Gms	Glass Bottle

### Available Packages

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com