



34321

Aluminium Oxide Activated (Basic) (Alumina)

Part A

CAS : 1344-28-1

Molecular Formula : Al2O3

Molecular Weight : 101.96

Specifications

Appearance (Colour)	Almost white
Appearance (Form)	Fine powder
pH (10% aq. suspension)	9.5 - 10.5
Chloride (Cl)	max. 0.01%
Sulphate (SO4)	max. 0.01%
Iron (Fe)	max. 0.01%
Heavy Metals (Pb)	max. 0.01%
Activity Grade (Brockmann)	1
Water Soluble Matter	max. 0.2%
Adsorption capacity (o-nitroaniline)	min. 1mg/g
Particle size (60 -325 Mesh BSS)	min. 90%
Pore size	58°A

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
25 Kg	28182010 (GST 18%)
5 Kg	28182010 (GST 18%)
500 Gms	28182010 (GST 18%)
Type of Packing	
5 Kg	Plastic Bottle
500 Gms	Plastic Bottle

Available Packages

500 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srlchem.com