



34516

Isobutylboronic Acid extrapure, 98%

Part B

CAS : 84110-40-7

Molecular Formula : C₄H₁₁BO₂

Molecular Weight : 101.94

Specifications

Appearance (Colour)	White to off - white
Appearance (Form)	Crystalline powder
Assay	min. 98%
Melting Point	105- 108°C

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:2709 , IMCO Class No.:3 , Packing Group:III
HSN Code	
1 Gms	29319090 (GST 18%)
25 Gms	29319090 (GST 18%)
5 Gms	29319090 (GST 18%)
Type of Packing	
25 Gms	Glass Bottle

Available Packages

1 Gms

5 Gms

25 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com