



34537

**L-Arginine Ethyl Ester Dihydrochloride  
extrapure, 98%**

Part B

CAS : 36589-29-4

Molecular Formula : C8H18N4O2.2HCl

Molecular Weight : 275.20

**Specifications**

Appearance (Colour)	White to off-white
Appearance (Form)	Crystalline powder
Assay (AT)	min. 98%

**General Information**

Storage	2 to 8°C (Refrigerate)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
100 Gms	29379020 (GST 18%)
25 Gms	29379020 (GST 18%)
5 Gms	29379020 (GST 18%)
Type of Packing	
100 Gms	Glass Bottle
25 Gms	Glass Bottle

**Available Packages**

5 Gms
25 Gms
100 Gms

**Disclaimer**

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srlchem.com