



34570

## Potassium Tetrafluoroaluminate pure, 98%

Part A

CAS : 14484-69-6

Molecular Formula : KAIF<sub>4</sub>

Molecular Weight : 142.07

### Specifications

Appearance (Colour)	Greyish
Appearance (Form)	Crystalline compound
Assay	min. 98%
Potassium (K) content	27 - 36%
Fluoride (F) content	44 - 54%
Aluminium (Al) content	16 - 21%

### General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
5 Kg	28261200 (GST 18%)
500 Gms	28261200 (GST 18%)
Type of Packing	
500 Gms	Plastic Bottle

### Available Packages

500 Gms  
5 Kg

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com