



34641

Aluminium Potassium Sulphate Dodecahydrate extrapure, 99%

Part A

CAS : 7784-24-9

Molecular Formula : $AlK(SO_4)_2 \cdot 12H_2O$

Molecular Weight : 474.38

Specifications

Appearance (Colour)	White
Appearance (Form)	Crystalline powder / Crystals
Solubility (Turbidity) 10% aq. solution	Clear
Solubility (Colour) 10% aq. solution	Colourless
Assay (Complexometric)	min. 99%
Chloride (Cl)	max. 0.005%
Iron (Fe)	max. 0.01%
Heavy Metals (Pb)	max. 0.002%
Ammonium (NH ₄)	max. 0.05%

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
2.5 Kg	28333030 (GST 18%)
500 Gms	28333030 (GST 18%)
Type of Packing	
500 Gms	Plastic Bottle

Available Packages

500 Gms

2.5 Kg

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com