



35161

1,3-Dimethyl-2-Imidazolidinone (DMI) GC-HS,
99.9%

Part A

CAS : 80-73-9

Molecular Formula : C5H10N2O

Molecular Weight : 114.15

Specifications

Appearance (Clarity)	Clear
Appearance (Form)	Liquid
Assay (GC)	min. 99.9%
Residual solvents	
GC-HS Class 1	max. 1ppm
GC-HS Class 2	max. 10ppm
GC-HS Class 3	max. 50ppm
Density (g/ml) @ 20°C	1.055 - 1.057
Refractive Index (20°C)	1.471-1.472
Boiling Point	224°C
Water (KF)	max. 0.1%

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
100 ml	29333100 (GST 18%)
500 ml	29333100 (GST 18%)

Available Packages

100 ml
500 ml

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srlchem.com