



35228

Silica Cyano Bonded extrapure

Part B

Other Information

Application

Can be used both in normal and reversed-phase chromatography as its polarity marks the separation between the polar and non-polar phases.

General Information

Storage

2 to 8°C (Refrigerate)

HSN Code

10 Gms

29319090 (GST 18%)

5 Gms

29319090 (GST 18%)

Available Packages

5 Gms

10 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srlchem.com