



35270

## Dextran ex. Leuconostoc Sp. - Tech 110 (100-120)

Part B

CAS : 9004-54-0

Molecular Formula : (C<sub>6</sub>H<sub>10</sub>O<sub>5</sub>)<sub>n</sub>

Molecular Weight : 1,10,000

### Specifications

Appearance (Colour)	White to slightly off-white
Appearance (Form)	Powder
Solubility (Turbidity) 1% aq. solution	Clear
Solubility (Colour) 1% aq. solution	Colourless
Specific rotation [?]	+195° to +203° (C=2 , Water)
Loss on drying	max. 7%
Average Molecular weight	100 - 120 KDa

### General Information

Storage	2 to 8°C (Refrigerate)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
100 Gms	39139090 (GST 18%)
25 Gms	39139090 (GST 18%)
5 Gms	39139090 (GST 18%)
500 Gms	39139090 (GST 18%)
Type of Packing	
25 Gms	Glass Bottle
5 Gms	Glass Bottle

### Available Packages

5 Gms  
25 Gms  
100 Gms  
500 Gms

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com