



35443

Magnesium Acetate Tetrahydrate extrapure, 98%

Part A

CAS : 16674-78-5

Molecular Formula : C₄H₆MgO₄.4H₂O

Molecular Weight : 214.40

Specifications

Appearance (Colour)	White
Appearance (Form)	Crystalline powder
Solubility (Turbidity) 10% aq. solution	Clear
Solubility (Colour) 10% aq. solution	Colourless
Assay (EDTA Titration)	min. 98%
Chloride (Cl)	max. 0.005%
Sulphate (SO ₄)	max. 0.01%
Heavy Metals (Pb)	max. 0.002%

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
5 Kg	29152920 (GST 18%)
500 Gms	29152920 (GST 18%)
Type of Packing	
5 Kg	Plastic Jar
500 Gms	Plastic Bottle

Available Packages

500 Gms
5 Kg

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srlchem.com, website www.srlchem.com